

14 September 2023

Unique Id: **SET-17508-A**

STATEMENT OF COMPLIANCE

to European Council Directive 2011/65/EU (recast of Directive 2002/95/EC) and 2015/863/EU on Restriction of the use of certain Hazardous Substances in electrical and electronic equipment.

The RoHS Directives require that electrical and electronic equipment placed on the EU market do not contain hazardous substances in excess of prescribed maximum concentrations within homogenous material, as below:

<u>Hazardous Substance</u>	<u>Maximum Concentration Level</u>
<u>Lead (Pb)</u>	<u>0.1% by weight (1000 ppm)</u>
<u>Mercury (Hg)</u>	<u>0.1% by weight (1000 ppm)</u>
<u>Hexavalent chromium (Cr VI)</u>	<u>0.1% by weight (1000 ppm)</u>
<u>Cadmium (Cd)</u>	<u>0.01% by weight (100 ppm)</u>
<u>Poly Brominated Biphenyl (PBB)</u>	<u>0.1% by weight (1000 ppm)</u>
<u>Poly Brominated Diphenyl Ether (PBDE)</u>	<u>0.1% by weight (1000 ppm)</u>
<u>Bis(2-ethylhexyl) phthalate (DEHP)</u>	<u>0.1% by weight (1000 ppm)</u>
<u>Butyl benzyl phthalate (BBP)</u>	<u>0.1% by weight (1000 ppm)</u>
<u>Dibutyl phthalate (DBP)</u>	<u>0.1% by weight (1000 ppm)</u>
<u>Diisobutyl phthalate (DIBP)</u>	<u>0.1% by weight (1000 ppm)</u>

Note: Products may qualify for exemption to the above limits, as defined in Annex III of RoHS2 Directive 2011/65/EU. Application-specific exemptions have not been applied.

Dynamic Controls confirms, to the best of our knowledge, that the items specified below, as sub-assemblies of electronic and electrical equipment supplied for further integration by a medical device manufacturer, conform to applicable requirements of Directive 2011/65/EU, (recast of Directive 2002/95/EC) and Directive 2015/863/EU.

The medical device manufacturer is responsible for final compliance of their finished device with regard to substance restrictions when equipped with the above item/s. Directive 2011/65/EU requires that finished electrical & electronic products be CE marked - components must not bear the CE Mark. This statement represents current knowledge following reasonable inquiry and is based on information provided by our suppliers and selective testing by independent laboratories.

14 September 2023

Unique Id: **SET-17508-A**

Rhino 2 Scooter Controller

Item	Dynamic Controls Part Number	Description
Scooter Controller	DS90	RHINO2 SCOOTER CONTROLLER 90 AMP
	DS120	RHINO2 SCOOTER CONTROLLER 120 AMP
	DS160	RHINO2 SCOOTER CONTROLLER 160 AMP
	DS180	RHINO2 SCOOTER CONTROLLER 180 AMP
	DS70C	RHINO2 CLASSIC CONTROLLER 70 AMP
	DS110C	RHINO2 CLASSIC CONTROLLER 110 AMP
	DS160C	RHINO2 CLASSIC CONTROLLER 160 AMP

For and on behalf of Dynamic Controls



Steve Beckett
Quality Manager