

Dynamic Controls
39 Princess Street, Riccarton 8041.
PO Box 1866, Christchurch 8140.
New Zealand
P: +64 3 962 2519
W: dynamiccontrols.com

15 July 2025

DECLARATION OF CONFORMITY

To:
European Radio Equipment Directive, 2014/53/EU (RED)

LiNX Access Keys

Part Number	Description
DLX-HKEY01-A	LiNX Access Key Dealer
DLX-HKEY02-A	LiNX Access Key OEM

Dynamic Controls declares that the products specified above, as service and diagnostic tool for use with electrically powered wheelchairs (Class I medical device), conform to all applicable requirements of the Radio Equipment Directive, 2014/53/EU.

Compliance with the relevant clauses of Article 3 – Essential requirements of RED 2014/53/EU is demonstrated by use of an integrated Bluetooth module that is fully compliant with RED 2014/53/EU Essential Requirements for Health, Safety, and Radio according to the following standards:

Radio: EN 300 328 v2.1.1

Health: EN 62479:2010

Safety: EN 60950-1:2006+AMD1:2010+ A2:2013


EMC: EN 301 489-1 v2.2.0, EN 301 489-17 v3.2.0

And final product testing to meet EN 301 489-1 v2.2.3

Dynamic Controls is an ISO13485:2016 Registered Supplier. Certified by SGS on 15 September 2023, Certificate AU02/55954.

This declaration of conformity is issued under the sole responsibility of Dynamic Controls.

For and on behalf of Dynamic Controls



Beth Gordon
Engineering Delivery Manager
Dynamic Controls
Christchurch NEW ZEALAND
Email: beth.gordon@allient.com

